



## CERTIFICATE OF ANALYSIS

### IHC Reference Standard

#### PRODUCT INFORMATION

Catalogue No.	HD720	Date of Manufacture	Mar 2015
Serial No.	13958	Format	FFPE Section on Slide
Description	BRAF V600E		

#### PROPERTIES

Storage Conditions	+4°C
Use By	Stable for 6 months from date of manufacture Use within 3 months of opening

#### QUALITY CONTROL

TEST	TEST METHOD	PASS/FAIL
Qualification	Protein expression confirmed using biomarker specific IHC assay	Pass
	Negative and positive cell line cores assessed using Quantitative Digital Pathology*	Pass

\* + Core (intermediate) shows weak staining intensity, thus statistical significance could not be obtained by QDP. A qualitative visual inspection was performed to confirm the identity of this core.

#### CUSTOMER SUPPORT

Technical	technical@horizondiscovery.com
Orders	orders@horizondiscovery.com

FOR RESEARCH USE ONLY