



CERTIFICATE OF ANALYSIS

IHC Reference Standard

PRODUCT INFORMATION

Catalogue No.	HD720	Date of Manufacture	August 2015
Serial No.	16011	Format	FFPE Section on Slide
Description	B-Raf V600E		

PROPERTIES

Storage Conditions	+4°C
Use By	Stable for 6 months from date of manufacture Use within 3 months of opening

QUALITY CONTROL

TEST	TEST METHOD	PASS/FAIL
Qualification	Protein expression confirmed using biomarker specific IHC assay*	Pass

* QDP analysis is not currently available.

CUSTOMER SUPPORT

Technical	technical@horizondiscovery.com
Orders	orders@horizondiscovery.com

FOR RESEARCH USE ONLY